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BORT AktiVen® donning aid for medical compression stockings

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a donning aid for medical compression stockings.

Indications

Ailments or disabilities with considerable functional disorders of the finger, hand or arm joints, primarily the gripping function which make the correct putting on of compression stockings prescribed by a physician impossible without aids, e.g. in case of highly inflammatory or considerably degenerative joint diseases, paralysis, deforming diseases in the hand region, results of injury, amputations; malformations, special sock or stocking donning aids for compression stockings with grip extensions are indicated due to the impossibility to reach to the stockings with the arms, e.g. in cases of far-reaching spinal, hip and knee joint stiffness as well as obesity

Contraindications

No contraindications are known if used properly.


Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.


- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- do not make any changes to the product
- do not use it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use – this medical aid is intended for treating one patient

Instructions for use:


In order to ensure successful use of your donning aid from the beginning, we kindly request you to read these instructions for use attentively. Use the BORT AktiVen® donning aid for medical compression stockings together with a table and chair with four freestanding chair legs.


Put the BORT AktiVen® donning aid for medical compression stockings on the table. Whilst holding the shell with one hand, pull the frame to its greatest possible length with the other hand until the side locking mechanism engages.


Place the shell in an upright, vertical position and the frame at the rear of the table.


Grasp the compression stocking with both hands at the upper edge so that the heel points towards the body.

Then slip the compression stocking evenly over the shell until the stocking's heel area is around the upper edge of the shell. Recommendation: For easier and energy-saving placement of the stocking, use special BORT AktiVen® gloves.

5
Set the splint hook so that the distance between the hook and the read edge of the shell is approximately as long as your foot.

6
Put the shell down on the base plate and set the frame upright.

7
Sit down on the chair with sufficient space in the immediate vicinity and place the donning aid with the base plate on the floor. Anchor the hook around the leg of a chair, securing the donning aid against forwards movement.

8
Continue to hold the frame and insert the foot into the compression sock/stocking from the rear to the heel.

Press your foot further forwards against the resistance of the compression stocking.

As soon as the shell straightens up, press your foot down until your heel contacts the floor.

Press your heel onto the floor and pull the frame in your direction until the upper edge of the shell reaches your popliteal fossa.

Place the frame aside.

Pull the compression stocking to the desired position.

If required, the knitted fabric can still be ideally spread over the leg and any creases in the compression stocking smoothed out. It's best to use the special BORT AktiVen® gloves for this step.

VIDEO

Material composition

ABS, light metal

Cleaning information

To clean, wipe with a damp cloth. Do not use any chemical cleaning agents. Neither grease nor oil the articulated joints.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the notified body for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EC) NO. 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 12.2019

 Medical device |  Single patient – multiple use

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