BORT Fußheberorthese

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
Med

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
Illustrationen
Figures


1


2


3

4

5

Sprachen Languages

DE	Deutsch	Gebrauchsanweisung	04
EN	English	Instructions for use	09
FR	Français	Mode d'emploi	14
ES	Español	Instrucciones de uso	19
IT	Italiano	Istruzioni per l'uso	24
NL	Nederlands	Gebruiksaanwijzing	29
CS	Čeština	Návod k použití	34
ET	Eesti	Kasutusjuhend	38
PL	Polski	Instrukcja użytkowania	42
RO	Română	Instrucțiunile de utilizare	47

PDF: ga.bort.com

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Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

Intended purpose

This medical device is a tibialis anterior muscle brace consisting of elastic material and rigid textile components. Stabilisation elements and correction straps support the tibialis anterior muscle in case of partial paralysis.

Indications

Peroneal paralysis, especially for mild to medium paralysis of the tibialis anterior muscle and the flexor hallucis longus (degree of severity 2–4).

Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

Application risks/Important notes

This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.

- remove the medical device prior to radiological examinations
- if wearing at night is prescribed, avoid negative influence on the circulatory system
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used

- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the brace, please neither use any creams nor ointments on or around it as they can destroy the material

Fitting instructions

Preparation for putting the BORT tibialis anterior muscle brace on:

The tibialis anterior muscle brace can be worn directly on the skin and over a thin sock. It's best to put the brace on when sitting.

1. Open all the tibialis anterior muscle brace Velcro fasteners and pull the brace over your foot.
2. Grasp the finger loop attached to the reverse side of the tibialis anterior muscle brace on the upper material edge and pull the tibialis anterior muscle brace upwards over the heel.

Closing the Velcro fasteners and attaching the restraints:

Before closing the Velcro constraints, please observe that the tibialis anterior muscle brace is correctly positioned on the foot. The tibialis anterior muscle brace must be correctly placed on the heel. For an effective course of the restraint, the seams must be central on the instep and in the heel area.

Important: During the entire following attachment process for the tibialis anterior muscle brace, the foot and heel remain on the floor!

- 1 Place your foot on the floor. Place restraints A and B on the outer side and restraint C on the inner side of the foot on the floor.
- 2 Then close the Velcro fasteners on the lower leg in the following order:
 - Firstly, the lower of the three fasteners.
 - In the next step, the fastener on the upper material edge.
 - Finally the fastener positioned in the middle.

If necessary, adjust the tension of the fasteners on the lower leg part of the tibialis anterior muscle brace again in order to achieve a pleasant but stable hold of the tibialis anterior muscle brace on the lower leg.

3 Slide the thigh or the knee as far forwards as possible. A foot angle of 20–30° to the lower leg is ideal. Guide restraint A upwards vertically first, then under tension diagonally over the instep and above the inner ankle around the lower leg. Velcro the restraint to the hook strap attached on the reverse side.

4 Guide the non-elastic restraint B upwards on the outer side of the foot under tension. At the same time, raise the outer edge of the foot slightly with the restraint. Then, attach restraint B to the velour Velcro surface of the tibialis anterior muscle brace.

Important: Restraint B must run upwards vertically. It stabilises the ankle joint and is intended to hinder twisting outwards when walking.

5 Finally, close restraint C. First, guide restraint C vertically upwards under tension diagonally over the instep and above the lateral malleolus around the lower leg. Then, attach restraint C to the hook ribbon attached to the rear side on restraint A. Attach restraint C to the restraint A hook strap as congruently as possible in order to avoid damages to socks.

Removal

To remove the BORT tibialis anterior muscle brace, open the belts in the reverse order.

Tip: Close the Velcro fasteners in order to avoid early wear to the Velcro surface.

Material composition

Cotton (CO), polyamide (PA), elastodiene/latex (LA), viscose (CV), elastane (EL)
The sewn in textile label provides the precise material composition.

The product contains latex and can trigger allergic reactions.

Cleaning information

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer

 Do not iron  Do not clean chemically

Do not use fabric conditioner. Stretch back into shape and dry in the fresh air. Close the Velcro fastenings to avoid damage to other laundry items.

Guarantee

The legal regulations of the country in which you acquired the product apply to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: www.bort.com/md-eu-contact.

Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: www.bort.com/conformity

Status: 10.2025

 Medical device |  Single patient – multiple use

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