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## EN BORT Stabilo Knee Support with Hinge

Many thanks for placing your trust in a medical device from BORT GmbH. Please read the existing instructions for use carefully. If you have any questions, please refer to your physician or the specialist retailer from whom you purchased this medical device.

### Intended purpose

This medical device is a brace for supporting the physiological guidance of the knee joint with circumferential Velcro fastener restraints and side guidance splints with joints.

### Indications

Mild to moderate lateral ligament instability of the knee joint, gonarthrosis, arthritis.

### Contraindications

Thrombosis risk, extreme varicosis, peripheral arterial occlusive disease (PAOD), lymph drainage disorders, also unclear soft tissue swellings distal to the aid positioned, sensory and circulatory disorders of the region of the body treated, skin diseases in the part of the body treated, latex allergy.

### Application risks/Important notes

- This medical device is a prescribable product. Discuss the use and duration with your treating physician. The expert staff from whom you have received the medical device will select the appropriate size and instruct you regarding its use.
- remove the medical device prior to radiological examinations
- remove it during longer periods of rest, e.g. when sleeping
- in case of numbness, loosen the medical device and remove it if necessary
- in case of persistent complaints, consult the physician or specialist retailer
- use the medical device in accordance with therapeutic needs
- only use other products simultaneously after consultation with your physician
- do not make any changes to the product
- do not wear it on open wounds
- do not use in case of intolerance of one of the materials used
- no re-use. This medical aid is intended for treating one patient only
- whilst wearing the support, please neither use any creams nor ointments on or around it as they can destroy the material

### Fitting instructions for BORT Stabilo knee support with articulated joint (only sizes 1 to 6+):

#### Information for the technician only:

Before putting the support on, the joint splints must be adapted to the contour of the body by an expert. The support is equipped with a joint for the inside leg (medial) and a joint for the outer leg (lateral). For forming, remove the fastening straps from the joint and take the splints out of the support. Then, form the joint splints to the contour of the leg so that the joint splint contacts the surface and causes no pressure sores. After this, the splints are inserted back in the splint bags and the straps adapted, but not closed fully.

#### 1 Information for the patient and technician:

- Information regarding individual length setting of the Stabilo strap using the easyClick system.
- ① attachment module, ② adjustment module. Adapt the length of the Stabilo strap individually by moving the adjustment module ②.
- When you receive your support, the attachment module ① of the Stabilo strap is attached to the support in the adjustment module ②. Firstly, open the attachment module ① on the lower leg and thigh by pulling and then hanging it out. Then, slip the support over the leg so that the patella is placed centred in the integrated silicone ring. Before closing the easyClick system, the length of the adjustment module ② must be adapted as described. Then, close the straps on the lower leg, followed by the thigh, by hanging the free attachment module ① onto adjustment module ②.
- Ensure that the Stabilo strap is placed comfortably and is neither too loose nor too tight.
- The support should now sit comfortably on the leg and not cause pressure sores. Slight corrections to the fit of the splints can also be made without removing them. **Important:** If the splints cause pressure sores and erythema, please remove the support and have it checked by the service provider.

### Fitting instructions for BORT Stabilo knee support with articulated joint (children's size only):

Before putting the support on, the joint splints must be adapted to the contour of the body by an expert. The support is equipped with a joint for the inside leg (medial) and a joint for the outer leg (lateral). For forming, remove the fastening straps from the joint and take the splints out of the support. Then, form the joint splints to the contour of the leg so that the joint splint contacts the surface and causes no pressure sores. After this, the splints are inserted back in the splint bags and the straps adapted, but not closed fully.

Then, slip the support over the leg so that the patella is placed centred in the integrated silicone ring. Then, close the straps on the lower leg, followed by the thigh. The support should now sit comfortably on the leg and not cause pressure sores. Slight corrections to the fit of the splints can also be made without removing them. **Important:** if the splints cause pressure sores and erythema, please remove the support and have it checked by the service provider.

#### Removal for BORT Stablio knee support with articulated joint (only sizes 1 to 6+):

To remove the support, open the easyClick system and hang it out. Grasp the lower edge of the support on the side splints and pull it downwards over the foot.

#### Removal for BORT Stablio knee support with articulated joint (children's size only):

To remove the support, open the Velcro fasteners, grasp the lower edge of the support on the splints and pull it downwards over the foot.



#### Material composition

Polyamide (PA), elastodiene/latex (LA), elastane (EL), viscose (CV)  
The sewn in textile label provides the precise material composition.

 The product contains latex and can trigger allergic reactions.

#### Cleaning information

Before washing, remove the joint splints and close the Velcro fastener in order to avoid damaging other items of clothing.

 Delicate wash 30 °C  Do not bleach  Do not dry in a tumble dryer  Do not iron  Do not clean chemically  
Do not use fabric conditioner. Stretch back into shape and dry in the fresh air.

#### Guarantee

The legal regulations of the country in which you acquired the product apply

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to the purchased product. Please contact your specialist retailer if you suspect a warranty claim. Please clean the product before submitting a warranty claim. If the enclosed instructions for use have not been properly observed, the warranty may be impaired or cancelled. The warranty does not cover use of the product inappropriate for the indication, non-observance of application risks, instructions and unauthorised modifications to the product.

#### Useful life/Lifetime of the product

The lifetime of the medical device is determined by its natural wear and tear if treated appropriately and as recommended.

#### Duty of notification

If a serious worsening of the state of health occurs when using the medical device, you are obliged to notify the specialist dealer or us as the manufacturer and the Medicines & Healthcare products Regulatory Agency (MHRA).

You can find our contact information in these instructions for use. You can find the contact information for the national authority responsible for your country under the following link: [www.bort.com/md-eu-contact](http://www.bort.com/md-eu-contact).

#### Disposal

Upon the termination of use, the product must be disposed of in accordance with the corresponding local requirements.

#### Declaration of conformity

We confirm that this device conforms with the requirements of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL. You can find the current declaration of conformity under the following link: [www.bort.com/conformity](http://www.bort.com/conformity)

Status: 10.2025

 Medical device |  Single patient – multiple use

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