
X Choice

 **Juzo**

JuzoFlex

Manu Xtra

Handgelenkbandage · Wrist support

Gebrauchsanweisung · Instructions for use · Mode d'emploi · Istruzioni d'uso · Instruções de utilização · Instrucciones de uso · Gebruiksaanwijzing · Brugsanvisning · Bruksanvisning · Käyttöohje · Bruksanvisning · Οδηγίες χρήσης · Kullanım talimatları · Upute za uporabu · Instrukcja użycia · Инструкция по эксплуатации · Návod na používanie · Navodila za uporabo · تعليمات الاستخدام

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EN

Intended purpose

JuzoFlex Manu Xtra provides anatomical guidance for the affected wrist, thus relieving it during movement. The adjustable plastic rod and the textile strap provide mechanical stabilisation. Medical compressive active supports have a proprioceptive effect. They improve self-perception and can compensate for muscle imbalances. They alleviate pain and support the removal of fluids from the tissue thanks to their rheological effect.

How to put on the product correctly:

1. Open the velcro strap and slide your hand into the wrist support until the stabilising splint lies on the inner surface of your hand. Put your thumb through the opening provided.
2. Position the wrist support in such a way that the ring-shaped depression in the pressure pad encloses the wrist bone without pressure.
3. Now close the velcro strap with sufficient tightness so that you feel an additional stabilisation of your wrist.

The correct fit of the wrist support is assured when it encloses the wrist tightly but does not apply a too strong (painful) pressure to the wrist.

The wrist support can be readjusted at any time by repeating steps 2 and 3 above.

Caution: rough, unfilled finger nails or sharp-edged rings etc. may damage the knitted fabric.

Please note:

Only wear your support when prescribed by your physician. Medical aids should only be dispensed by suitably trained staff of a medical products supplier. The support can only be fully effective when it has been put on correctly and is being worn in accordance with treatment recommendations. The support is largely resistant to oils, ointments, skin moisture and other environmental influences. The durability of the material may be affected by certain combinations of the above. We therefore recommend that it be examined regularly by the medical products supplier. Should the product become damaged, please consult your specialist supplier. Do not repair the support yourself – this could compromise its quality and medical efficacy. For hygienic reasons, the support is not suitable for treating more than one patient. The sewn-in textile label is important for the identification and traceability of the product. We therefore ask you not to remove the label under any circumstances.

Washing and care instructions

Please note the care instructions on the textile label. Please remove the reinforcing stay from the support. Remove the strap and close it to protect the hook and loop areas from dirt and damage. Wash the supports separately the first time (colour may run). Your support should be washed every day or between uses on a gentle or delicate wash cycle (30 °C). Tip: a laundry net provides additional protection for the knitted fabric. Please do not use fabric softener! When hand washing, rinse well and do not wring. We recommend the gentle Juzo special detergent. You can reduce drying time by placing the support on a thick towel after washing, tightly rolling the towel up and firmly pressing out excess moisture. Then hang the support up loosely. Do not leave the product in the towel and do not dry it on a radiator or in the sun. Juzo supports are not suitable for dry cleaning. The support is not suitable for tumble drying. The removed stabilising splint can be cleaned with a damp cloth. Then slide it back into the pocket provided for it when the support is completely dry (with the bent spoon in front – the open side of the hollow section faces away from the hand surface). To finish, fasten the strap to the support.

Material composition

For precise details, please see the textile label that is sewn into your product.

Storage information and usage duration

The medical retailer or doctor prescribing the support will recommend exactly how long it should be used for.

If the product is being worn and washed every day, its maximum usage duration will be 6 months. If the product is used on an irregular basis, the wearing duration will be longer, up to a maximum of 12 months. Based on individual use, product wear may differ. We therefore recommend that the product is checked by a medical specialist in the event of suspected changes. Store the product in a dry place and protect it from sunlight. Supports have a maximum shelf life of 42 months. The usage duration of the medical product is printed on the box with an hour-glass symbol. Body measurements may change based on the patient's condition and individual circumstances. We therefore recommend that body measurements are checked regularly by a medical specialist.

Indications

Distortions (sprains / strains), tendosynovitis (tendovaginitis), arthrosis (degeneration of the wrist), arthritis (inflammation of the wrist), post-traumatic irritation of the wrist.

Contraindications

If you have any of the following conditions, you should only wear the support after consulting with your doctor: Lymph drainage disorders and non-specific soft tissue swelling outside the area of application, sensory disturbances and impaired circulation in parts of the body to be treated, skin conditions, skin conditions, skin irritation and injuries in the area to be treated; open wounds in the treatment area must be covered with sterile material. Liability for failure to comply with these contraindications cannot be accepted by the Julius Zorn GmbH Group.

Side effects

There are no known side effects when using this product correctly. However, if you notice any adverse changes (e.g. irritation of the skin) during the prescribed use, please consult your doctor or medical retailer immediately. If you know that you are allergic to one or more of the product ingredients, please consult your doctor before using the product. If your symptoms worsen during the time you are wearing the product, please take off the product and consult a medical specialist immediately. The manufacturer does not accept any liability for damage/injuries resulting from improper use or misuse. In the case of complaints related to this product, such as damage to the knitted fabric or problems with the fit, please contact the medical retailer directly. Only serious incidents that may lead to a major deterioration of the patient's medical condition or to death should be reported to the manufacturer and the relevant authority in the Member State. Serious incidents are defined in Article 2 (65) Regulation (EU) 2017/745 (MDR).

Disposal

As a rule, the product is to be destroyed via the conventional disposal route after the end of use. Please observe the applicable local/national regulations.

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**Julius Zorn GmbH** · Juliusplatz 1 · 86551 Aichach · Germany · info@juzo.de · juzo.com
**Berro AG** · Austrasse 7 · 9055 Bühler · Suisse
UK responsible person: Juzo UK Ltd