
50R49 Dorso Carezza Posture

DE	Gebruuchsanweisung	3
EN	Instructions for Use	7
FR	Instructions d'utilisation	12
IT	Istruzioni per l'uso	16
ES	Instrucciones de uso	21
NL	Gebruiksaanwijzing	26
PL	Instrukcja użytkowania	31
HU	Használati utasítás	35
CS	Návod k použití	40
SK	Návod na používanie	44
TR	Kullanma talimatı	49
RU	Руководство по применению	53
CN	使用说明书	58

	Size/ Größe	Waist circumference/ Taillenumfang*	
		cm	inch
	XXS	60 – 70	23.6 – 27.6
	XS	70 – 80	27.6 – 31.5
	S	80 – 90	31.5 – 35.4
	M	90 – 100	35.4 – 39.4

*approx. 2 cm/0.8 inch above the iliac crest / ca. 2 cm oberhalb des Beckenkamms

Material	Fabric/ Gewebe: PA, EL, Lycra®
	Corset Splints/ Korsettstäbe: POM

1 Verwendung

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1 Application

1.1 Intended Use

The orthosis is intended **exclusively** for orthotic fittings of the trunk and **exclusively** for contact with intact skin.

1.2 Effects

The orthosis activates muscular stabilisation of the thoracic spine and straightens it actively. It promotes sensorimotor functions, supports pain relief, and relieves the thoracic spine.

1.3 Indications

- Poor body posture (hunchback without fused vertebrae, Scheuermann's disease)
- Back pain (thoracic spine)

Indications must be determined by the physician.

1.4 Contraindications

1.4.1 Absolute Contraindications

None.

1.4.2 Relative contraindications

Consult the physician if any of the following occur: skin diseases/injuries, inflammatory changes, raised scars with swelling, reddening, and warming in the area where the orthosis is fitted, tumours, lymphatic congestion – including unclear swelling of soft tissue distal to the device, and dysaesthesia in the legs.

1.5 Further usage restrictions

The orthosis is designed for use on **one** patient. Re-use is not permissible. The orthosis must be used according to the indications and following the instructions of a physician.

2 Safety

2.1 Explanation of warning symbols

 WARNING	Warnings regarding possible risks of severe accident or injury.
 CAUTION	Warnings regarding possible risks of accident or injury.
 NOTICE	Warnings regarding possible technical damage.

2.2 General Safety Instructions

 CAUTION
The support / orthosis may be over-tightened
Risk of local pressure and constriction of blood vessels and nerves.
► Verify the correct positioning and fit of the support / orthosis.

CAUTION

Use on more than one patient

Skin irritation, formation of eczema or infections

- Only use the support / orthosis on one patient.

NOTICE

Incorrect environmental conditions

The support / orthosis is not flameproof.

- Do not expose the support / orthosis to open flame or other heat sources.

NOTICE

The lifespan of the materials may be reduced

- Avoid contact with substances containing oils and acids, salves and lotions.

INFORMATION

Improper changes to the support / orthosis are not permitted and will eliminate product liability.

2.3 Information to be Provided to the Patient

CAUTION

Direct skin contact with the support/orthosis can cause skin irritation

- Inform the patient that he/she should not wear the support/orthosis directly against the skin.

NOTICE

Damage due to wear and tear

- ▶ Instruct the patient to inspect the support / orthosis for proper functionality, wear and damage prior to each use.
- ▶ Advise the patient that the support / orthosis must not be used if a closure or other part of the support / orthosis shows signs of wear.

INFORMATION

The patient must be instructed to see a physician immediately if any exceptional changes are noted (e. g. worsening of the complaint).

3 Handling

INFORMATION

The daily duration of use and usage period are generally determined by the physician. Initial fitting and application of the support/orthosis must be performed by trained specialists. The patient must be instructed on how to handle and care for the support/orthosis.

3.1 Selecting the size

- 1) Measure the waist circumference approx. 2 cm above the iliac crest.
- 2) Determine the orthosis size (see size table).

3.2 Fitting and Application

- 1) **Application** (see Fig. 1): Place the orthosis over the shoulders like a backpack.
- 2) **Hook and loop material:** Route the ends of the lower strap to the front (see Fig. 2) and secure the hook and loop closure (see Fig. 3). Route the ends of the shoulder strap to the front around the chest and secure the hook and loop closure below the chest (see Fig. 4)

- 3) **Final check:** The correct fit and proper positioning of the orthosis must be checked during handover to the patient.
- 4) **Instruct the patient:** Explain to the patient how to don the orthosis and ensure that it is seated correctly. Practice donning and doffing of the orthosis with the patient.

3.3 Cleaning

Clean the orthosis regularly.

- 1) Close the hook and loop closures.
- 2) **Recommendation:** Use a laundry bag or net.
- 3) Wash the orthosis in 40 °C warm water with standard mild detergent. Do not use fabric softener.
- 4) Allow to air dry. Do not expose to direct heat (e.g. from the sun, an oven or radiator).

4 Disposal

Dispose of the product in accordance with national regulations.

5 Legal information

5.1 Liability

The manufacturer's warranty applies only if the device has been used under the conditions and for the purposes described. The manufacturer recommends that the device be used and maintained according to the instructions for use.

5.2 CE Conformity

This device meets the requirements of the 93 / 42 / EEC guidelines for medical devices. This device has been classified as a class I device according to the classification criteria outlined in appendix IX of the guidelines. The declaration of conformity was therefore created by Otto Bock with sole responsibility according to appendix VII of the guidelines.

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Otto Bock HealthCare GmbH