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INSTRUCTIONS FOR USE EN

Dear patient,

we ask you to carefully follow the enclosed instructions for use. If you have any questions, please contact the doctor looking after you, your nearest specialist store or us directly.

INTENDED USE

EPI-HiT® is an elbow brace for relieving the muscle origins at the ulnar and / or radial epicondyle.

INDICATIONS

- Conservative / postoperative
- Epicondylitis humeri radialis/ulnaris (tennis elbow, golfer's elbow)

CONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swelling, redness) of the body areas to be supplied
- Impaired circulation or lymphatic soft tissue swelling
- Neurogenic disorders of the sensory system and skin trophism in the body area to be treated (sensory disorders with or without skin damage)

SIDE EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indica-

tions and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE

- The exact positioning of the EPI-HiT® brace is approx. 3-4 cm below the inflammatory changed tendon tissue or directly on the pathologically hypertonic (tense) muscle. Make sure that the large friction pad is applied to the affected side.
- Adjust the brace via the latch mechanism by lifting the tab slightly with the fingernail **1** so that its ends lie on the muscle bellies, see illustration for correct fit. Adjust the pad accordingly.
- Now close the Velcro fastener.
- The fastener of the brace is continuously adjustable and the brace should be adjusted with it so that a firm fit is achieved without uncomfortable pressure.
- The correct pressure is achieved when the brace can no longer be felt after a few minutes when the arm is relaxed, but pressure is immediately felt under the pad when the hand is used.
- If the brace is too tight, loosen it slightly or move it down a few millimetres. Also move the brace downwards if it pinches when the elbow is bent in flexion.
- The brace should be worn all day, please take it off at night.
- The wrist support should surround the wrist positively. With the stabilising band, you can now adjust the degree of support as needed.

NOTES FOR THE ORTHOPAEDIC TECHNICIAN

- A closure strap that is too long can be shortened in the Velcro area at one of the welds. The smaller pad should lie in the overlap area under the larger pad. The Velcro point can be moved if necessary. The pad can be shortened in the narrow area at one of the thin embossed strips if necessary.

ACCESSORIES / SPARE PARTS

Optionally available from your specialist store:

- Replacement pad for EPI-HiT® ([REF](#) 37505)
- Wrist support ([REF](#) 01125)

⚠ IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- If a significant deterioration in health occurs during the use of the medical device, we ask you to report this serious incident to your specialist dealer, your doctor or us as the manufacturer, as well as to the responsible authority. You can find our contact information in these instructions for use.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.

- Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.
- SPORLASTIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

CARE

We recommend washing the EPI-HiT® brace and wrist support gently in lukewarm water (30 °C) with mild detergent and air drying. Do not dry the product in the tumble dryer, on the heater or in the microwave. It is not necessary to remove the pads for washing. Close the Velcro fasteners when washing to keep them functional for longer and to avoid damage to other garments.

⚡ TRANSPORT AND STORAGE CONDITIONS

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product

at a normal temperature and humidity.

QUALITY MANAGEMENT SYSTEM

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

An online version of the instructions for use can be found on our website at www.sporlastic.de.

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