REF ART.-NR. 07642
MANU-X®

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INSTRUCTIONS FOR USE EN

Dear Patient,

We ask you to carefully follow the enclosed instructions for use. If you have any questions, please contact the doctor looking after you, your nearest specialist store or us directly.

INTENDED USE

MANU-X® is a wrist brace for immobilisation of the wrist in one direction or plane of movement with free movement of the fingers.

INDICATIONS

- Conservative/postoperative
- Arthritis of the wrist
- Severe sprains
- Early functional postoperative follow-up after the first dressing change
- Overuse arthropathies
- Interim treatment after plaster casts
- Conservative treatment of initial stages of aseptic bone necrosis of the carpus
- Drop hand

KONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swelling, redness) of the body areas to be supplied.
- Impaired circulation or lymphatic soft tissue swelling.
- Neurogenic disorders of the sensory system and skin trophism in the body area to be treated (sensory disorders with or without skin damage).

SIDE EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE

- Slip your hand into the MANU-X® wrist brace **1A**.
- Place the Velcro strap around the wrist and fasten it with the desired pull **2A, 3A**.
- If necessary, you can detach the Velcro from the micro Velcro again, readjust it and fasten it.
- MANU-X® has the correct fit when the volar reinforcement splint, which is preformed in the neutral position, fits positively against the palm, wrist and forearm **4A**.
- MANU-X® should have a firm, stabilising fit on your hand without pressing uncomfortably.

NOTES FOR THE ORTHOPAEDIC TECHNICIAN

- Right (R) Left (L) marking on the wrist brace **1B**.
- By turning the wrist brace 180 degrees, MANU-X® can be used for the left or right hand. Depending on the side of application, the enclosed reinforcement splint on the inside of the hand must be inserted accordingly.
- The R/L marking with the stylised hands shows you the position for the application of MANU-X®.
- When applying the MANU-X® for the right hand (marked R), the hand on the marking points with the thumb to the left **1B**.
- When applying for the left hand (marked L), the hand on the marking points with the thumb to the right.
- When inserting the reinforcement splint, make sure that the short curved side always points towards the fingertips **2B**.
- Insert the splint accordingly into the slot made on the inside of the brace **3B**. Make sure that it is then enclosed on both sides **4B**.
- The reinforcement splint on the back of the hand is already integrated. Its position does not have to be changed when putting on the brace. If necessary, the splint

can be removed through a slot and individually adjusted.

IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- In the event of serious incidents when using this aid, immediately contact your doctor, your health care supply retailer or contact SPORLASTIC directly at info@sporlastic.de to report the incident.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.
- Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.
- SPORLASTIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

CARE

We recommend washing MANU-X® gently in lukewarm water (30°C) with mild detergent and air drying. Do not dry the product in the tumble dryer, on the heater or in the microwave. For washing, please remove the reinforcement splints from the brace. Close the Velcro fasteners before washing. These remain functional longer and damage to other items of clothing can be avoided. Please do not use fabric softener when washing the MANU-X®.

TRANSPORT AND STORAGE CONDITIONS

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

QUALITY MANAGEMENT SYSTEM

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

- An online version of the instructions for use can be found on our website at www.sporlastic.de.

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Kõik SPORLASTIC GmbH tooted läbivad Saksaama kvali- tedijuhitumise süsteemi nõuetele vastava ülevaatuse.