	<p>DE Trocken aufbewahren vor Feuchtigkeit schützen</p> <p>EN Keep dry protect against moisture</p> <p>FR Conserver au sec Protéger à l'abri de l'humidité</p> <p>IT Conservare in un luogo asciutto lontano dall'umidità</p> <p>PL Przechowywać w suchym miejscu i chronić przed wilgocią</p> <p>DA Opbevares tørt beskyttes mod fugt</p>
	<p>DE Vor Sonnenlicht schützen</p> <p>EN Protect from sunlight</p> <p>FR Protéger des rayons du soleil</p> <p>IT Tenere al riparo dalla luce solare</p> <p>PL Chronić przed promieniowaniem słonecznym</p> <p>DA Beskyttes mod sollys</p>
	<p>DE Einzelner Patient, mehrfach anwendbar</p> <p>EN Single patient, multiple use</p> <p>FR Patient unique, application polyvalente</p> <p>IT Singolo paziente, uso multiplo</p> <p>PL Jeden pacjent, wielokrotne użycie</p> <p>DA En patient, kan anvendes flere gange</p>
	<p>DE Webseite zur Patienteninformation</p> <p>EN Patient information website</p> <p>FR Site Internet d'informations pour les patients</p> <p>IT Sito web con informazioni per i pazienti</p> <p>PL Strona internetowa z informacjami dla pacjentów</p> <p>DA Hjemmeside til patientoplysninger</p>
	<p>DE Nach Gebrauch wegwerfen</p> <p>EN Discard after use</p> <p>FR Usage unique. Jeter après utilisation</p> <p>IT Smaltire dopo l'uso</p> <p>PL Wyrzucić po użyciu</p> <p>DA Bortskaf efter brug</p>
	<p>DE Verwendbar bis</p> <p>EN Usable until</p> <p>FR Date d'expiration</p> <p>IT Utilizzabile fino a</p> <p>PL Data ważności</p> <p>DA Kan bruges indtil</p>
	<p>DE Hersteller des Produktes</p> <p>EN Manufacturer of the product</p> <p>FR Fabricant du produit</p> <p>IT Produttore dell'articolo</p> <p>PL Producent produktu</p> <p>DA Producent</p>
MD	<p>DE Medizinprodukt</p> <p>EN Medical device</p> <p>FR Dispositif médical</p> <p>IT Prodotto medico</p> <p>PL Urządzenie medyczne</p> <p>DA Medicinsk produkt</p>
REF	<p>DE Artikelnummer</p> <p>EN Article number</p> <p>FR Référence</p> <p>IT Numero articolo</p> <p>PL Numer artykułu</p> <p>DA Varenummer</p>
LOT	<p>DE Fertigungslosnummer, Charge</p> <p>EN Production batch number, batch</p> <p>FR Numéro de lot de fabrication</p> <p>IT Numero di lotto di produzione, partita</p> <p>PL Numer partii produkcyjnej, szarża</p> <p>DA Batchnummer, batch</p>
PZN	<p>DE Pharmazentralnummer</p> <p>EN Pharmaceutical registration number</p> <p>FR Numéro central de la pharmacie</p> <p>IT PZN (Germania)</p> <p>PL Numer centralny farmacji</p> <p>DA Farmacentralnummer</p>
SPORLASTIC (GmbH)	<p>DE SPORLASTIC (GmbH) ist als Kurzform der Firmenbezeichnung SPORLASTIC GmbH Medizinische Produkte zu verstehen</p> <p>EN SPORLASTIC (GmbH) is to be understood as a short form of the company name SPORLASTIC GmbH Medizinische Produkte</p> <p>FR SPORLASTIC (GmbH) est une forme abrégée de la dénomination sociale SPORLASTIC GmbH Medizinische Produkte</p> <p>IT SPORLASTIC (GmbH) è da intendersi come forma abbreviata della ragione sociale SPORLASTIC GmbH Medizinische Produkte</p> <p>PL SPORLASTIC (GmbH) należy rozumieć jako skróconą formę nazwy firmy SPORLASTIC GmbH Medizinische Produkte</p> <p>DA SPORLASTIC (GmbH) skal forstås som en forkortelse af firmanavnet SPORLASTIC GmbH Medizinische Produkte</p>

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INSTRUCTIONS FOR USE EN

INTENDED USE

NEURODYN® CARBON / NEURODYN® CARBON NRG is a foot drop brace to support a physiological gait pattern in the case of existing foot drop weakness.

INDICATIONS

[REF] 07853, NEURODYN® CARBON:

- Conservative / postoperative
- Mild poliomyelitis
- Flaccid and (moderately severe) spastic foot drop paresis (e.g. after stroke, multiple sclerosis, peroneal nerve palsy, other neurological/muscular diseases)
- Impairment of walking in the event of neuromuscular damage to the foot lifting muscles (flaccid paralysis in peroneal nerve palsy, CNS damage, nerve root damage to the lumbar spine)

[REF] 07854, NEURODYN® CARBON NRG:

- Conservative / postoperative
- Mild poliomyelitis
- Flaccid and (moderately severe) spastic foot drop paresis (e.g. after stroke, multiple sclerosis, peroneal nerve palsy, other neurological/muscular diseases)
- Forefoot amputation (partial)
- To correct the pointed foot position/lifting of the foot, active return of the foot to the individual functional range of motion and maintenance of mobility
- Partial paralysis of the foot lifting muscles (and high mobility of the insured person or the insured person) due to e.g. peroneal nerve palsy, lumbar root damage, stroke

CONTRA-INDICATIONS

- Allergic, inflammatory or injury-related skin changes (e.g. swellings, redness) of the body areas to be cared for
- Impaired circulation or lymphatic soft tissue swelling
- Neurogenic disorders of the sensory system and skin trophism in the body area to be cared for (sensory disorders with and without skin damage)
- Strong supination position
- Lymphatic drainage disorders
- Circulatory disorders, PNP
- Severe spastic foot drop paresis

SIDE-EFFECTS

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

INTENDED PATIENT TARGET GROUP

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

PUTTING IT ON AND USE

- Remove the original sole of your shoe **1**.
- Open both straps at the Velcro end that sticks to the lower leg shell. It has an orange ring **2**.
- Place NEURODYN® CARBON / NEURODYN® CARBON NRG in the shoe and make sure that the foot plate fits snugly with the glued-on insole **3**.
- Slip into the shoe in a sitting position. After you have put on the shoe, position your foot forwards until your lower leg is in contact with the frontal support **4**.
- Now grasp the open Velcro end of the upper fastener and guide it around the lower leg. Fasten the Velcro end with the attached deflection loop on the strap area that lies over the frontal attachment **5**.
- Repeat this process with the lower strap **6**.
- If necessary, you can now use the Velcro ends with the black ring to fine-tune the belt tension **7, 8**.
- In the area of the lower leg frame, the course of the straps is controlled by Velcro points. If the inclination of the strap does not fit the lower leg perfectly, you can lift the strap off the shell and reposition it.

Selection of the shoe:

- The selection of the right shoe is crucial for the function of the NEURODYN® CARBON / NEURODYN® CARBON NRG.
- Discuss the correct footwear with your technician. The effective heel height of the shoe should be approx. 10 mm.
- Closed shoes with a removable footbed are best suited. This provides sufficient volume for the brace. It is important that the heel is enclosed.
- There is no body weight limit for the brace.
- There are also no restrictions for everyday activities.
- However, the brace is not approved for sporting activities.

NOTES FOR THE ORTHOPAEDIC TECHNICIAN

- NEURODYN® CARBON NRG / NEURODYN® CARBON can neither be ground nor thermoformed.
- For fitting in the shoe, the brace size is determined according to the shoe size.
- The carbon sole of the orthosis is always smaller than the specified shoe size. The insole can be cut to fit the shape of the original insole **1A**.
- The position of the orthosis in the shoe is determined and secured by gluing on the insole. The double-sided adhesive strip on the cover sole can be used for gluing.
- Use the enclosed insole to determine the pressure-free position of the brace in the shoe. Pressure points can only be eliminated by pushing the orthosis away from the pressure point. To do this, the insole must be removed and reattached in the new position. The orthosis can then not slide back towards the pressure point due to the form-fitting insole in the shoe.
- Localised pressure points along the edge of the shinbone can also be made hollow by gluing the enclosed padding strips to the back of the shell pad.
- Ensure that the orthosis is not slanted in the shoe if the heel is too high and that it is not slanted backwards if the heel is too low. Check the heel spring of the shoe and adjust if necessary.
- To check, it is best to place the shoe with the brace in it next to the contralateral leg. The front edge of the frontal support should now be at the same angle as the tibia of the contralateral leg when standing upright **2A, 3A**.
- If necessary, readjust the heel height of the orthosis using the levelling cut-outs included in the scope of delivery **4A**.
- The bar pad included in the scope of delivery ensures that the required minimum distance is maintained between the bar and the back of the foot during fitting.
- You should always check the success of the fitting by carrying out a dynamic fitting. Please allow the user to walk over uneven terrain and inclined surfaces. The fitting result should also be checked while seated.
- Instead of the insole, we recommend using a customised orthopaedic soft foam insole.

ACCESSORIES

Optionally available from your specialist retailer:

- Hygiene set consisting of shell pad, bar pad and two fastening straps with fastening rivet for NEURODYN® CARBON / NEURODYN® CARBON NRG (REF 17853)

IMPORTANT INFORMATION

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- If a significant deterioration in health occurs during the use of the medical device, we ask you to report this serious incident to your specialist dealer, your doctor or us as the manufacturer, as well as to the responsible authority. You can find our contact information in these instructions for use.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Please check NEURODYN® CARBON / NEURODYN CARBON® NRG at regular intervals. If you notice any wear, abrasion, etc., please visit your medical supply retailer immediately. We generally recommend that you visit your medical supply retailer approximately every 6 months to check the brace.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.
- To drive a motor vehicle/machine safely, the driver must be fit to drive. The responsibility for this lies with the vehicle/machine operator. In the event of an accident, you may be required to prove that you were able to drive the vehicle/machine safely. If necessary, you can have your suitability checked and confirmed by an authorised body.

 Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.

 SPORLASTIC GmbH is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

CARE

We recommend cleaning the splint and shin shell of the NEURODYN® CARBON / NEURODYN® CARBON NRG with a damp cloth and gently washing the padding in lukewarm water (30°C) with mild detergent and air drying. The padding of the shin guard can be unclipped for cleaning. Do not dry the product in a tumble dryer, on a radiator or in the microwave.

TRANSPORT AND STORAGE CONDITIONS

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

QUALITY MANAGEMENT SYSTEM

All products from SPORLASTIC GmbH are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

An online version of the instructions for use can be found on our website at www.sporlastic.de.

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