de Webseite zur Patienteninformation
en Patient information website
fr Site Internet d'informations pour les patients
it Sito web con informazioni per i pazienti
pl Strona internetowa z informacjami dla pacjentów
da Hjemmeside til patientoplysninger
sl Spletno mesto z informacijami za bolnike

de Nach Gebrauch wegwerfen
en Discard after use
fr Usage unique. Jeter apres utilisation
it Smaltire dopo l'uso
pl Wyrzucić po użyciu
da Bortskaf efter brug
sl Po uporabi zavržite

de Verwendbar bis
en Usable until
fr Date d'expiration
it Utilizzabile fino a
pl Data ważności
da Kan bruges indtil
sl Uporabno do

de Hersteller des Produktes
en Manufacturer of the product
fr Fabricant du produit
it Produttore dell'articolo
pl Producent produktu
da Producent
sl Proizvajalec izdelka

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Instructions for use

Dear patient,

We ask you to carefully follow the enclosed instructions for use. If you have any questions, please contact the doctor looking after you, your nearest specialist store or us directly.

Intended use

VERTEBRADYN® X.FORCE is a lumbar spine brace with mobilisation function for relief by deloading the lumbar spine from L1 to L5.

Indications

All indications where relief through decompression followed by functional mobilisation of the lumbar spine is necessary, such as:

- Severe lumbar sciatica with muscle deficits in cases of significant disc protrusion/prolapse, without absolute indication for surgery
- Severe radicular, pseudo-radicular lumbar syndrome (resistant to conservative therapy)
- Grade II spondylolisthesis with recurrent lumbar sciatica
- Severe lumbar deformity with facet syndrome/arthritis
- Condition after disc surgery without residual paresis with cleared disc
- Spinal canal decompression, one day, without deformity, disc surgery one day
- Spinal canal stenosis with short walking distance without significant paresis

All indications where relief of the posterior vertebral edge is necessary, such as:

- Minor post-traumatic vertebral fractures
- Pathological fractures, e.g. fractures with risk of plasmacytoma
- Inflammation
- Degeneration
- Damage in the lumbar spine
- For lumbar spine decompression with the option of adjustment for early functional expansion of the range of motion

- Lumbar instability
- Hyperlordosis pain with significant severity
- Scheuermann's disease of the lumbar spine in adults with symptoms
- Osteoporosis: isolated lumbar spine fracture without deformity with pain symptoms

Contra-indications

- In general, the treating physician should be consulted regarding the indication and wearing mode of an orthopaedic aid, both in general and in the presence of the following medical conditions in particular:
- Allergic, inflammatory or injury-related skin changes (e.g. swelling, redness) in the areas of the body to be treated
- Impaired circulation or lymphatic soft tissue swelling
- Neurogenic sensory and skin trophic disorders in the body area to be treated (sensory disorders with and without skin damage)

Side effects

If used appropriately and correctly, there are no general side effects, intolerances with long-term effects or allergic reactions known to date.

Intended patient target group

The target group is all patients, taking into account indications and contraindications. Instruction and allocation of the correct size is carried out by medical professionals.

Putting it on and use

- We recommend wearing a thin, close-fitting T-shirt under the product.
- Open the pull elements and fix them loosely in front **1**.
- Open the belly button closure. Use the hand straps to put on the brace and guide it around your body **2**. Position the back section centrally over the spine. The orange lower edge of the brace should end above the sacrum. The frame should extend from the lower ribs to the waistline.
- Close the abdominal fastener while applying tension so that the brace provides comfortable pressure and the abdominal plate fits centrally **3**.
- Tighten the two lower pull elements evenly and secure them to the abdominal fastener **4**. Repeat this with the upper pull elements **5**.
- Make sure that the front lower edge is flush and parallel to the thigh when sitting and does not cause any pressure in the groin **6**.
- To remove the brace, open the straps and loosely secure the pull elements **7**.
- Open the abdominal fastener and remove the brace.
- VERTEBRADYN® X.FORCE requires individual and qualified adjustment to the patient by appropriate specialist personnel. This is the only way to ensure performance and wearing comfort.

Notes for the orthopaedic technician

- The height of the brace frame can be adjusted. To do this, remove the brace frame from the frame padding **A**.
- Use a 2.5 mm Allen key to loosen the corresponding screws and adjust the height **B**.
- The wings of the lower frame section should be positioned between the iliac crest and the greater trochanter. The wings of the upper frame section should enclose the lower rib **C**.
- Once the height has been adjusted correctly, secure the screws in the appropriate positions on the brace frame **D**.

- The back frame is made of hard aluminium and can be shaped by hand or with a bending iron .
- Secure the brace frame with frame padding in the centre of the brace .
- The spring steel rods can be shaped individually.
- For increasing mobility, the brace can be gradually dismantled in consultation with the attending physician. To do this, the frame or the abdominal plate can be removed from the brace.
- With the removal of the abdominal plate and frame, the VERTEBRADYN® X.FORCE becomes a lumbar support brace.

Accessories

Optionally available from your specialist retailer:

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Important information

- Essentially, the indication for and duration of wearing of any orthopaedic device in general, as well as the presence of any of the conditions below in particular, should be discussed with the doctor who is in charge of your treatment.
- Check the functioning of the product with your doctor or specialist retailer.
- The product is designed to treat one patient.
- Do not wear the product directly against open wounds.
- In order for the product to have a long service life and function, it may not be worn in conjunction with oily/greasy or acidic products, ointments or lotions.
- The product can come into contact with splash water.
- Please always close any Velcro fasteners present properly. Damage caused by improper closing of the Velcro fasteners does not constitute a reason for complaint.
- If a significant deterioration in health occurs during the use of the medical device, we ask you to report this serious incident to your specialist dealer, your doctor or us as the manufacturer, as well as to the responsible authority. You can find our contact information in these instructions for use.
- Any improper modifications to the product and/or improper use of the above-mentioned product exclude any product liability on the part of the manufacturer.
- Product may contain metallic objects. We recommend not putting them in contact with electromagnetic interference.
- Possible health-related mutual risks or other disadvantages with certain treatments that may arise in association with the use of the product must be discussed with the doctor in charge of treatment.

 Please return the packaging and the product to the local recycling collection system. Please comply with local regulations.

 SPORLASTIC GmbH & Co. KG is affiliated with a returns system for packaging and therefore fully complies with the packaging legislation in the Federal Republic of Germany.

Care

We recommend washing VERTEBRADYN® X.FORCE gently in lukewarm water (30 °C) with mild detergent and air drying. Do not dry the product in the tumble dryer, on the heating or in the microwave. Remove the frame before washing. The cover and the abdominal plate can be washed together with the frame. Close the Velcro fasteners during washing to ensure that they remain functional for longer and to prevent damage to other items of clothing. Washing restores the tension of the moulded knit, which

diminishes over time. We therefore recommend washing the brace more frequently.

Transport and storage conditions

Please ensure that the product is stored in a dry place and is protected against moisture and sunlight. Store the product at a normal temperature and humidity.

Quality management system

All products from SPORLASTIC GmbH & Co. KG are subject to product inspection as part of our quality management system. If you still have any complaints about our product, we kindly ask that you contact your specialist retailer.

Verwendetes Material

Polyamide, polyvinyl chloride, polyester, polyurethane, elastane.

 An online version of the instructions for use can be found on our website at www.sporlastic.de.

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