
MANU-CAST®

ORGANIC FX P LONG

SPORLASTIC

REF 07859

MANU-CAST® ORGANIC FX P LONG

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for use

Dear Patient,

This orthopaedic device is classified as a medical device, so we request that you carefully read these Instructions for use. If you have any questions, please contact your attending physician, your nearest specialist dealer or us directly.

Intended use

The MANU-CAST® ORGANIC FX P LONG is a wrist-hand-thumb brace for immobilisation and restraint of the wrist in two directions of movement with thumb fixation and free movement of the fingers.

Application/indications

All indications in which restraint or immobilisation is necessary following damage to the wrist/middle hand area and thumb, with free movement of the fingers (pre-operative/post-operative/post-traumatic/conservative), e.g.:

- Stable distal radius/ulna fractures
- Stable scaphoid fractures
- Stabilisation and immobilisation of the first metacarpophalangeal joint (MCP)
- Stabilisation and immobilisation of the carpometacarpal joints
- After removal of plaster casts
- Carpal tunnel syndrome/CTS
- Inflammations
- Degeneration
- Rheumatic conditions
- Distortions
- Disc injury
- Acute wrist injuries
- Osteoarthritis of the wrist
- Radial nerve palsy
- Kienbock's disease
- Navicular pseudo-arthritis

Contraindications

- Conditions involving restricted arterial or lymphatic circulation in the affected area
- Neurogenic sensitivity and trophic disorders of the skin in the affected area

Risks/side effects

At present, SPORLASTIC is not associated with any side effects, provided it is used properly and applied correctly. Signs of pressure and damage to blood vessels or nerves may occur nevertheless in connection with fixed devices. Please discuss the following with your attending physician before using the support:

- Illnesses in the area of application
- Sensory disturbances, circulatory disorders
- Lymphatic drainage disorders

Intended users and target patient group

Patients after assignment and guidance by medical professionals, with due consideration of the indications and contraindications. Healthcare professionals, including carers and relatives who look after the patient and who have received instruction from healthcare professionals.

Application and use

- Open the Velcro fasteners of the MANU-CAST® ORGANIC FX P LONG and pull the brace over the hand.
- The anatomically pre-shaped mesh frames are positioned on the inside and outside of the hand.
- Position the upper (dorsal) mesh frame on the outside of the hand so that it is centred.
- The upper bar must lie between the thumb and index finger.
- Close the Velcro fastener on the back of the hand with a measured pull ❶
- Next, fasten the wrist strap. This is wrapped circularly around the wrist and the brace ❷
- Now close the two forearm straps from the front (hand) to the back (forearm) ❸

- If necessary, retighten the circular wrist strap.
- Close the circular Velcro fastener on the thumb ❹
- The MANU-CAST® ORGANIC FX P LONG is correctly tightened when both anatomically pre-shaped mesh frames fit snugly on the palm and outer surface of the hand as well as on the thumb, wrist and forearm ❺
- To remove the MANU-CAST® ORGANIC FX P LONG, open the Velcro fasteners and pull off the brace with the other hand.

Reuse ❸

MANU-CAST® ORGANIC FX P LONG is intended exclusively for the individual treatment of a single patient.

Important notes ❶

- As a rule, the medical use and function of an orthopaedic device, as well as how long it is worn, should always be discussed with the attending physician or specialist retailer.
- Kindly refer also to the symbols on the packaging and these Instructions for use.
- In order to ensure a long service life, the product must not be worn in conjunction with fatty and acidic products, ointments or lotions.
- Please do not use any aggressive detergents or disinfectants.
- Do not wear the product on open wounds.
- If included, always ensure that all Velcro fasteners are properly closed. Damage caused by improper closing of the Velcro fasteners shall not constitute grounds for complaint.
- Consult with your attending physician or specialist dealer if you notice a deterioration in your state of health while using this device.

Disposal ❷

Please comply with local regulations when disposing of the packaging and product after use.

Care and maintenance

Gently wash the MANU-CAST® ORGANIC FX P LONG in lukewarm water (max. 40 °C) with mild detergent and then leave to dry without direct exposure to sunlight. Do not dry the product in a tumble dryer, on the radiator or in the microwave. The orange plastic clips can be opened to remove the pad from the brace for cleaning. Wash the product separately or close any Velcro fasteners first.

The product is essentially maintenance-free with the correct handling and care.

Transport and storage conditions ❹

Store the product at room temperature and away from moisture and direct sunlight.

Complaints

Kindly contact your specialist dealer if you have any complaints about our product. If you return the product to us, please clean it beforehand.

Notification obligations

Please report any physical damage due to a product defect to SPORLASTIC and your national authority for the supervision of medical devices. You will find our SPORLASTIC contact details at the end of these Instructions for use.

Disclaimer

Please use the medical device only as discussed with your attending physician or specialist dealer. No liability is accepted for self-medication. Unprofessional modifications to the product and/or improper use of the product shall void any manufacturer's product liability.

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❸ An online version of these Instructions for use is available on our website at www.sporlastic.de.

Instructions for the orthopaedic technician

- MANU-CAST® ORGANIC FX P LONG can be individually adjusted and readjusted as required.
- To remould, immerse the appropriate point of the mesh frame (ORGANIC FRAME) in hot water at 70-80 °C. Thermoplastic remoulding of the heated area is then possible.
- Take note: Do not use a hot air blower for remoulding, as heating will be uneven!
- The fasteners of the MANU-CAST® ORGANIC FX P LONG are infinitely adjustable on both sides.

Accessories and spare parts

The following spare parts are enclosed with the product or can be obtained as options from your specialist dealer:

- MANU-CAST® ORGANIC FX safety strap (REF 37459)
- Replacement thumb strap (REF 27159)

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